#### **Gmail Account**

thesisvitalsignsmapua@gmail.com

thesisvitalsignsmapua01

Note: Wag galawin yung Google Sheets sa Gmail Account, ayun ginamit na database. Pwede iexport, wag lang iedit.

#### **000 WebHost Account**

thesisvitalsignsmapua@gmail.com

thesisvitalsignsmapua01

**Website Link:** : (Hosted in 000 webhostapp free hosting, not available between 2 AM - 3 AM for free version)

https://vitalsignsmonitormapua.000webhostapp.com

### **Program Files Link:**

Website (For PC view only):


- https://github.com/romnegrillo/THESIS\_VITALIOTWEBB

Device Programs (In RPi and Arduino):


- https://github.com/romnegrillo/THESIS\_VITALIOT

## Manual

- First register sa website.
- Mga registered users lang makakagamit sa device.
- May guest mode pero di masasave sa website.
- Follow niyo na lang instruction sa device.








### **Website Manual**


Pag admin, kaya niya maaccess lahat ng data ng user.


- View user vital signs.
- View user symptoms
- Search specific vital signs for specific date.
- View/edit user personal data.
- Delete user.

Username: admin

Password: admin








## Pag user:

- View symptoms
- View vital signs



# Starting the device

- Pag may keyboard mouse kayo, rekta run niyo na lang program. Nasa folder na THESIS\_VITALIOT/Programs/vitalsigns\_main.py
- Check niyo muna kung may internet kayo, bago irun program kasi magsstuck lang yung sa isang windows since aattempt attempt niyo magconnect.
- Open niyo cmd sa RPi tas ping lang kayo kahit anong website.